

# Containment Test Data Task Specific

Containment Performance Target (CPT) = 10ng/m<sup>3</sup>

Test		No. Test Runs	Outcome (task based)		
Ref:	Device/Activity:		Personal	Static	Overall
02	Sample removal through continuous liner	3	3 of 3 results <10% CPT	9 of 9 results <10% CPT	PASS
03	Toxin linker dispensing and transfer	3	3 of 3 results <10% CPT	15 of 15 results <10% CPT	PASS
04	Power failure	3	3 of 3 results <10% CPT	15 of 15 results <10% CPT	PASS
05	Induced containment breach	3	3 of 3 results <10% CPT	15 of 15 results <10% CPT	PASS
06	Flexible isolator contaminated enclosure disposal	1	2 of 2 results <10% CPT	3 of 3 results <10% CPT	PASS (indicative)